

**DUAL
TECHNOLOGY
WELLNESS
SYSTEM**

**ELECTRO
FREQUENCY
+
HYPERLIGHT**

Fullerene C60 Optics

5 Red Hyperlight Modes

36 Gold-Coated Needle Grid

19 Microcurrent
Intensity Levels

Frellux

GET WELL · STAY WELL



zepter

A photograph of a woman's back and shoulder. Two red dotted circles are overlaid on the skin, one on the upper shoulder and one on the lower back, indicating treatment areas. The woman's hands are resting on her shoulder and back.

Frellux

GET WELL

STAY WELL

FRELLUX™ represents an **evolution in personal care**, uniquely combining **Electrofrequency (EF) Microcurrent** and **Red Hyperlight Frequency (LF)** within a **single system** for **effortless, fast and effective relief** from **tension, aches, and everyday discomfort** while supporting **natural tissue recovery** processes and **skin health**.

Anytime, anywhere.

Frellux™ features three selectable modes:

- 1. Electrofrequency Microcurrent Stimulation**
19 adjustable intensity levels for personalized comfort and depth

- 2. Red Hyperlight Frequency Stimulation**
5 dedicated programs designed for different applications and sensitivities

- 3. Electrofrequency + Red Hyperlight Frequency (Synergy Mode)**
Combined application for a more comprehensive and enhanced effect



TWO FREQUENCIES.
TWO MECHANISMS OF ACTION.
ONE DEVICE.

Red Hyperlight Frequency penetrates the tissue to optimize microcirculation and cellular activity. By enhancing oxygen and nutrient flow, it effectively harmonizes the biological response in stressed or overworked areas.

Electrofrequency Microcurrent acts deeper within the musculoskeletal structure, delivering gentle signals to stimulate muscle fibers and release built-up tension. It supports natural neuromuscular function and accelerates the body's recovery in fatigued areas.

Non-invasive. Effortless. Effective.



36 gold-coated micro-needle grid



ELECTRO-FREQUENCY MICROCURRENT STIMULATION

Frellux™ delivers microcurrent therapy through a **36 gold-coated micro-needle grid**, providing wide, evenly dispersed stimulation that ensures consistent delivery and reduces the need for precise point placement. The electrical stimulation interacts with the body's natural bioelectrical signaling, supporting communication between cells, encouraging local circulation, and helping relax tight or overworked muscle fibers—contributing to natural recovery and everyday physical comfort.

With **19 adjustable electrostimulation energy levels**, each session can be customized from ultra-gentle to deeply soothing, allowing the intensity to be matched precisely to comfort, sensitivity, and personal preference.



BENEFITS OF ELECTRO-FREQUENCY MICROCURRENT STIMULATION

Microcurrent stimulation supports local circulation, improving oxygen and nutrient delivery to tissues helping the body release tension faster and restore normal function in overworked areas.

- ⊕ **Less muscle tension and tightness**
- ⊕ **Reduced stiffness, especially after sitting or inactivity**
- ⊕ **Relief from everyday aches and discomfort**
- ⊕ **Easier, more comfortable movement**
- ⊕ **Less tired or heavy feeling in the body**
- ⊕ **Faster recovery after physical effort**




a quantum mechanical light transformer

HYPERLIGHT FULLERENE C₆₀

Any light passing through the hyperlight optics (lens) is transformed into a uniquely structured hyper-harmonized light complex, called quantum hyperlight.

C60 molecules, rotating at extremely high speeds, reorganize the light passing through them into a highly structured and symmetrical form. This structure reflects the symmetry associated with many healthy biomolecules.

Through the resonance principle of biomimicry, this structured Hyperlight interacts more effectively at the level of biomolecules, helping support the natural rebalancing of disturbed or stressed systems—promoting recovery, relieving tension, and helping maintain healthy tissue and skin function.





RED HYPERLIGHT FREQUENCY STIMULATION

Frellux™ utilizes a 640-660 nm Red Hyperlight wavelength, designed to interact with the skin and superficial tissues through photobiomodulation. This light supports natural cellular activity and healthy microcirculation, helping improve oxygen and nutrient delivery at the surface level.

As a result, it helps calm sensitive or irritated areas, support tissue comfort, and promote recovery after stress or strain. It also contributes to improved skin tone, smoother texture, and a more balanced, healthy appearance.

LF Mode offers five distinct programs from a Steady Continuous Hyperlight that supports circulation and helps promote deep tissue comfort and relaxation to High-Frequency Pulsation Hyperlight up to 8,000 Hz, designed for gentle stimulation of sensitive or inflamed areas without generating thermal stress.

Frellux

GET WELL * STAY WELL



5 INDIVIDUALLY SELECTABLE HYPERLIGHT PROGRAMS

PROGRAM	FREQUENCY	CORE CHARACTERISTICS	TYPICAL SUPPORT
1.	Continuous	Strong, steady illumination	Circulation and stiffness relief
2.	10 Hz	Slow rhythmic pulses	Muscle relaxation & comfort
3.	600 Hz	Revitalizing pulses	Circulation and collagen activity
4.	3000 Hz	Non-heating pulses	Sensitive or reactive tissue
5.	8000 Hz	Ultra-fine pulses	Surface comfort and delicate skin

BENEFITS

OF RED HYPERLIGHT

- ⊕ Softens stiffness in muscles
- ⊕ Eases localized musculoskeletal discomfort, supporting better mobility
- ⊕ Soothes irritation-prone or highly sensitive areas
- ⊕ Calms and soothes inflamed or reactive areas
- ⊕ Supports natural tissue repair processes
- ⊕ Enhances skin tone and natural radiance
- ⊕ Supports skin firmness and elasticity
- ⊕ Helps restore the skin after environmental or physical stress

Frellux™ is an ideal handheld, on-the-go solution for both physical comfort and cosmetic restoration.



SAFE & RESPONSIBLE USE



For optimal results, we recommend using the pre-set treatment time of 10 minutes, after which the device switches off automatically. Shorter sessions can be performed by switching the device off at any time.

During use, apply slow, circular movements to avoid localized overheating during longer sessions and to assure even energy distribution when treating larger areas.

Avoid Electro-Frequency (EF) Microcurrent use:

- On open wounds, inflamed, irritated, or broken skin.
- If you have a pacemaker or any implanted electrical or electronic medical device.
- During pregnancy, unless advised by a healthcare professional.
- Directly on the head, the front of the neck (thyroid area), or along the spine.
- In cases of active bleeding, bleeding disorders, or a known tendency to hemorrhage.
- In cases of thrombosis or known blood clotting disorders.

Avoid Red Hyperlight Frequency (LF) use if you have:

- Acute febrile conditions (fever or active infection).
- Known photosensitivity or light-sensitive conditions
- Epilepsy or a history of seizures

If you have a medical condition, are undergoing treatment, or are unsure whether this device is appropriate for you, consult a qualified healthcare professional before use.

Always refer to the instructions for use provided with the device for complete guidance on safe operation.





1. Electrofrequency Intensity levels: 1-19

2. Red Hyperlight program 1-5

3. Battery level

TECHNICAL INFORMATION

Frellux	AP-001
MADE IN	POLAND
WEIGHT	250G
CHARGING TIME	3 HOURS
IMPULSE INTENSITY ADJUSTMENT	YES - FROM 1 TO 19
APPLICATION	MAX 20 CONSECUTIVE MINUTES
TOTAL OPERATING TIME AFTER BATTERY CHARGE	8 CYCLES 2X10MIN. WITH A 3-MINUTE BREAK
OUTPUT VOLTAGE	5.0 V DC
CHARGER	USB-C 5.0V DC

zepter®

Frellux



HOME ART & SALES SERVICES AG
Sihleggstrasse 23, CH-8832 Wollerau, Switzerland,
Telefon: +41 43 888 28 00, Email: info@zepter.com.
Copyright HOME ART & SALES SERVICES AG - ZEPTER
INTERNATIONAL 2026 PMD-ME-159-27-EN

SHOP NOW
zepter.com

